



# Custom AFO Collection SELECTION GUIDE

## Ankle Stabilizers

Indications: Ankle instability, minimal hindfoot and midfoot control



Model:	SURE-01	SURE-02	SURE-03	SURE-08
Ankle control	Minimal	Moderate	Dorsiflexion assist	Selective ROM
Hind foot control	Minimal	Minimal	Minimal	Minimal
Midfoot control	Minimal	Minimal	Minimal	Minimal

## Ankle Gauntlets

Indications: Ankle instability, moderate hind foot and midfoot control



Model	SURE-19	SURE-16	SURE-18	SURE-14	SURE-13
Ankle control	Moderate	Selective ROM	Maximal	Maximal	Maximal
Hind foot control	Moderate	Moderate	Maximal	Maximal	Maximal
Midfoot control	Moderate	Moderate	Moderate	Moderate	Moderate

## Traditional AFOs

Indications: Ankle instability and hindfoot, midfoot, forefoot and knee control



Model	SURE-10	SURE-11	SURE-12	SURE-09	SURE-07
Ankle control	Minimal	Plantarflexion stop	Dorsiflexion assist	Maximal	Maximal
Hind foot control	Moderate	Moderate	Moderate	Maximal	Maximal
Midfoot control	Moderate	Moderate	Moderate	Moderate	Maximal
Forefoot control	Moderate	Moderate	Moderate	Moderate	Maximal
Knee	Minimal	Moderate	Moderate	Maximal	Maximal



# Custom AFO Collection SELECTION GUIDE

Condition	MILD	MODERATE	SEVERE
Ankle instability	SURE-01	SURE-02	SURE-14,18,13
Pes planus	SURE-01	SURE-02	SURE-14,18,13
Pes cavus	SURE-01	SURE-02	SURE-14,18,13
Peroneal nerve palsy	SURE-10	SURE-03,08,11,12,16	SURE-10,11
Achilles tendonitis	SURE-01,02		
Sinus tarsi syndrome	SURE-01	SURE-02	SURE-08,14,18
Drop foot	SURE-10	SURE-03,08,11	SURE-11,09
PTTD	SURE-01	SURE-02	SURE-08,16,12
Post Fracture	SURE-01,02	SURE 08,16	SURE-08,14,18,12
Tarsal coalition	SURE-01,02,08	SURE-14,16	SURE-08,09
Ulcers			SURE-07
ROM limiter		SURE 08,16	SURE-12
Charcot degeneration	SURE-01,02,10	SURE-14,18,11,12	SURE-13,09,07
Triplanar ankle foot control			SURE-08,09
Knee hyperextension	SURE-10	SURE-11,12	SURE-09
Excessive plantarflexion	SURE-03	SURE-08,16,18,14,13	SURE-11,09
Swing clearance	SURE-10	SURE-03,08,16	SURE-11,12,09
Poor Balance		SURE-19	



### STEP ONE: Indicate foot/ankle condition

Flexible

Fixed / Rigid



### STEP TWO: Include (5) key measurements

ML at Met heads

Finished brace height

ML at ankles

Circumference at finished brace height

Finished footplate length



### STEP THREE: Provide mold of foot and ankle

Cast model should represent final brace position (corrected foot position)

\* Simply, complete our order form. Select from a variety of models and options. All three steps must be met to qualify.  
\* 90 warranty on workmanship.